

## PCR and allele-specific restriction analysis (ASRA) of *pfdhfr* codons 50, 51 and 59

### PRIMARY PCR

5 ul of slide/filter-paper extract in total 25 ul PCR

Sense primer FR519-A 5'**GCGCGCTAATAACTACACATTTA**3'

Antisense primer FR519-B 5'**CCCGGGCTCTTATATTTCAATTT**3'

Product size: 147 bp

#### PCR PROGRAM (ASRA-1)

Primary Denaturing: 95°C, 5 min

Denaturing: 92°C, 30 sec

Annealing: 45°C, 30 sec

Extension: 65°C, 45 sec

Cycles: 45

Final Extension: 72°C, 15 min

### SECONDARY PCR FOR 51 AND 59 CODONS

1-2 ul of primary reaction in total 25 ul secondary PCR

Sense primer FR51-D

5'**CTAGGAAATAAAGGAGTATTACCATGGAAATGGA**3'

Antisense primer FR59-D

5'**ATTTTTTCATATTTTGATTCATTCACATATGTTGTAAGTGTAC**3'

Product size: 113 bp

#### PCR PROGRAM (ASRA-2)

Primary Denaturing: 95°C, 5 min

Denaturing: 92°C, 30 sec

Annealing: 45°C, 30 sec

Extension: 65°C, 30 sec

Cycles: 15-25

Final Extension: 72°C, 15 min

### SECONDARY PCR FOR 50 CODON

Sense primer FR519-A **GCGCGCTAATAACTACACATTTA**

Antisense primer FR59-D

**ATTTTTTCATATTTTGATTCATTCACATATGTTGTAAGTGTAC**

Product size: 128 bp

#### PCR PROGRAM (ASRA-2)

Primary Denaturing: 95,5 min  
Denaturing: 92, 30 sec  
Annealing: 45, 30 sec  
Extension: 65, 30 sec  
Cycles: 15-25  
Final Extension: 72, 15 min

#### RESTRICTION DIGEST

Digests for each were performed separately. 5-8 ul of PCR was digested with 1 unit enzyme overnight at 37°C in a total 20 ul reaction using New England Biolabs #4 buffer. Products were examined by agarose electrophoresis on 2% NuSeive (FMC, Rockville, MD), with an "uncut" (no enzyme) digest run alongside.

#### 51 codon

Enzyme: EcoRI (New England Biolabs, Beverly, MA)

Codon cleaved: wild (N)

Cleaved product sizes: 35 and 78bp

#### 59 codon

Enzyme: BsrGI (New England Biolabs, Beverly, MA)

Codon cleaved: wild (R )

Cleaved product sizes: 43 and 65bp

#### 50 codon

Enzyme: Tai I (New England Biolabs, Beverly, MA)

Genotype cleaved: mutant (Arg) and not wild (Cys)

Cleaved product sizes: approx. 58 and 70 bp

## PCR and allele-specific restriction analysis (ASRA) of *pfdhfr* codons 108 and 164

### PRIMARY PCR

Sense primer FR100-A **GGGGGGCAGTTACAACATATGTGA**

Antisense primer FR100-B **GGGGGCACATTCATATGTACTATTT**

Product size: 414bp

#### PCR PROGRAM (ASRA-1)

Primary Denaturing: 95,5 min

Denaturing: 92, 30 sec

Annealing: 45, 30 sec

Extension: 65, 45 sec

Cycles: 45

Final Extension: 72, 15 min

### SECONDARY PCR

Sense primer FR108-D

**CTAATTCTAAAAAATTACAAAATGT**

Antisense primer FR164-D3

**TTTCTTTTCTAAAAATCTTGATAAACAACGGAACCTCTTA**

Product size: 254 bp

#### PCR PROGRAM (FR100-2)

Primary Denaturing: 95,5 min

Denaturing: 92, 30 sec

Annealing: 42, 30 sec

Extension: 65, 45 sec

Cycles: 15-25

Final Extension: 72, 15 min

### RESTRICTION DIGEST

Digests for each are performed separately

NOTE: for Psi I, use 0.5X NEB#4 buffer and 1x BSA

Analyze products on 2% NuSeive gel

#### 108 codon

Enzyme: Alu I (New England Biolabs, Beverly, MA)

Codon cleaved: wild (Ser) and not mutants (Asn or Thr)

Approx. cleaved product sizes: 46 + 210 bp

Note: BsrI will cut 108 Asn; ScrFI will cut 108 Thr.

#### 164 codon

Enzyme: Psi I(New England Biolabs, Beverly, MA)

Codon cleaved: wild (Ile) and not mutant (Leu)

Approx. cleaved product sizes: 42 and 214bp

## PCR and allele-specific restriction analysis (ASRA) of *pfdhfr* codon 16

### PRIMARY PCR

Sense primer AMP-1 TTTATATTTTCTCCTTTTTA

Antisense primer FR519-B 5'CCCGGGCTCTTATATTTCAATTT3'

Product size: 255 bp

#### PCR PROGRAM (ASRA-1)

Primary Denaturing: 95,5 min

Denaturing: 92, 30 sec

Annealing: 45, 30 sec

Extension: 65, 45 sec

Cycles: 45

Final Extension: 72, 15 min

### SECONDARY PCR

Sense primer SP1 ATGATGGAACAAGTCTGCGAC

Antisense primer FR59-D

5'ATTTTTCATATTTTGATTCATTCACATATGTTGTAAGTGTAC3'

Product size: 220 bp

#### PCR PROGRAM (ASRA-2)

Primary Denaturing: 95°C, 5 min

Denaturing: 92°C, 30 sec

Annealing: 45°C, 30 sec

Extension: 65°C, 30 sec

Cycles: 15-25

Final Extension: 72°C, 15 min

### RESTRICTION DIGEST

5-8 ul of PCR was digested with 1 unit enzyme for 5 hours in a total 20 ul reaction using New England Biolabs #2 or 4 buffer. Analyze products on 2% NuSeive gel with an "uncut" (no enzyme) digest run alongside.

### 16 codon

Enzyme: Mwo I (New England Biolabs, Beverly, MA)

Codon cleaved: wild (Ala) and not mutant (Val)

Cleaved product sizes: 174 and 46 bp